

SEPTEMBER
2023

High throughput sequencing

Functional genomic data analysis: transcriptomics

Stéphane Le Crom
stephane.le_crom@sorbonne-universite.fr



First generation sequencing methods


Sanger sequencing by synthesis

Method discovered in 1977 by Frédéric Sanger (nobel price 1980).

DNA polymerisation using a **complementary primer**. Elongation using **thermostable DNA polymerase** (PCR).

Addition of 4 **deoxynucleotides** (dATP, dCTP, dGTP, dTTP) and low concentrations of one of four **dideoxynucleotides** (ddATP, ddCTP, ddGTP ou ddTTP).

These ddNTP once incorporated in the newly synthesized DNA strand, block elongation. Synthesis termination is done by a statistical manner on each possible positions.



Sequence reading


We get a **mix of DNA fragments** terminating at each position of the sequence.

These fragments are then separated on a DNA **polyacrylamide gel electrophoresis**.

Detection of synthesized fragments is done by the **incorporation of labelling beacon** in the DNA.

At the origin this label was radioactive, attached either on the primer or on the dideoxynucleotide.

Around **1 kb of DNA** by run during **2 days**.
One read by **sample**.





Capillary sequencers

First version in the 90's thanks to the **modification of the radioactive label by a fluorescent one**.

Using **glass capillary of few micron diameters**, on 30 to 50 cm long.

The four nucleotides migrate in the same tube thanks to **four different fluorescent dyes**.

300 kb of DNA by run during **3 hours**. Several hundred sample at a time.



ABI 3730xl DNA Analyzer

96 parallel capillary (up to 50 cm) array.
768 samples, 690 kb DNA, 3 hours run.


At the Broad Institute (Cambridge, Massachusetts) **126 devices** were able to sequence **1 human genome** in **12 days**.

HOW IT WORKS

Sequencing lanes create a long way back light. Back then, applied fluorescent dye attached DNA sequences. The graft, around bands ending of the lanes, it's also get "to" to each locally automated. Fluorescence-labeled and passed, for gets to band, and can the machine identify. In this, all lanes a good get would, some get "back get" sample get, so have time, so processing full would that genome delivery, so have get this.

Today, Applied Biosystems' sequencing system allows lanes to be run fully, a fluorescence-based system. Each surface in the separation column and gel, and the tube cap, are run get of them to parallel, generating some per base ends. So more exact, sample loading, the machine simply reads the capillary array, into a sample plate, into it, gets full automatic on the base and to read the sequencing reactions and process the sequence data, and a full single pump into a machine based on this.

So, presented, as the sequencing is, it's full automatic, the same in this. The final feature, in technology, from, mostly used in sequence time to call each component, through it, requires a better genome of one of this.



From *The Scientist*

Second generation sequencing: high throughput sequencing

The first technologies on the market

Goal: to obtain a huge number of short reads

since 1996



Applied Biosystems
ABI 3730XL
1 Mb / day



January '07 GS20
June '07 GS FLX

Roche / 454
Genome Sequencer FLX
100 Mb / run

January '07



Illumina / Solexa
Genome Analyzer
2,000 Mb (2 Gb) / run

November '07



Applied Biosystems
SOLiD
3,000 Mb (3 Gb) / run

Illumina Genome Analyzer

January 2007

```
From: Clive Brown <clive.Brown@solexa.com>  
Date: Sun, 20 Feb 2005 16:34:46 +0100  
To: Nick McCooke <Nick.McCooke@solexa.com>, Tony Smith <Tony  
Swerdlow <Harold.Swerdlow@solexa.com>, John Milton <JM.Milton  
<Kevin.Hall@solexa.com>, Colin Barnes <Colin.Barnes@solexa.c  
<Vincent.Smith@solexa.com>, Klaus Maisinger <Klaus.Maisinger  
Conversation: WE'VE DONE IT !!!!  
Subject: WE'VE DONE IT !!!!
```

Tony Cox, Peta and I now agree - having looked at all of the PhiX174 data.

We have re-sequenced our first genome !!!!!!!



DNA library preparation


Random DNA fragmentation and size selection.

Ligation of adaptors.

DNA denaturation.

Hybridization of fragments onto the “flowcell” surface.

Solid phase bridge PCR.




Reversible terminator sequencing

The four **reversible terminators** are added simultaneously.

Laser scanning of the flowcell surface.

Release of the blocking terminator.

Sequencing cycles are repeated one base at a time.




<http://www.illumina.com/>

Sequence analysis


Scanning at each position for all sequences (reads) in parallel.

Most of the errors (99%) are sequencing errors (misincorporation).



<http://www.illumina.com/>

Illumina sequencing by synthesis



<https://www.youtube.com/watch?v=fCd6B5HRaZ8>

Specifications of the latest Illumina sequencers




	MiniSeq	MiSeq	NextSeq 550	NextSeq 2000	NovaSeq
Run Time	24 hours	56 hours	29 hours	2 days	44 hours
Read length (bp)	2x 150	2x 300	2x 150	2x 150	2x 150
Read number	50 10 ⁶	50 10 ⁶	800 10 ⁶	1 10 ⁹	10 10 ⁹
Ouput	7.5 Gb	15 Gb	120 Gb	300 Gb	3,000 Gb
Throughput	7 Gb/day	7 Gb/day	100 Gb/day	150 Gb/day	1,500 Gb/day

The third generation

Real time sequencing

Real time sequencing on **single molecule** thanks to **RNA polymerase immobilisation** in wells.

Each **base incorporation** is measure in real time with a CCD camera under the bottom of the plate.



Sequel specifications

8,000,000 reads / SMRT cell;

From 1 to 52 kb (avg 20 kb);

450 Gb by SMRT cell;

Run duration = 30 hours.



Nanopore technology


Single molecule detection system by passing single strand nucleic acids through a **nanometric pore**.

Base to base analysis in real time using **electric properties** of the nanopore.

DNA size sequencing of **kilobases**.

No limitation on the acid nucleic type to be detected (**DNA or RNA** even amino acids) including modifications (epigenetics).

No amplification.



Greenwood (2013) *Popular Science*


First described by David Deamer in 1989

Sunday June 25, 1989

Driving back from Eugene → Belknap Lodge, had an idea on how to sequence DNA directly:

Main concept:
DNA will be driven through a small channel, either by $\Delta\psi$ or ΔpH . The channel will be carrying a current, driven by $\Delta\psi$. As each base passes through, a change in the current will occur. Because the bases are of different size, the current change will be proportional, thereby providing an indication of which base it is.

Details:
The thickness of the membrane must be very thin, perhaps a polymerized bilayer. The channel must be of the dimensions of DNA in cross section, approx. 1-2 nm. Porin? Complement? Alamethicin? The ion flux might be protonic.




Oxford nanopore technologies



The MiniON flow cell

1 flow cell = 1 membrane
with 512 nanopores.

Single molecule sequencing
up to 4 Mb during up to 72 hours.



@NanoporeConf


Field genomics: portable sequencer tracks infectious disease outbreaks



Nick Loman using a MinION to sequence the Zika virus in Brazil

Applications

They cover a lot of previous existing techniques



Kahvejian et al. (2008) *Nat. Biotech.*

De novo sequencing

Quicker and cheaper sequencing than Sanger.



But **reads are smaller**.

Combination of different methods allow to obtain **better quality sequencing drafts**.

=> Combining 454 and Illumina.

Low error rate and **homogenous coverage** due to no cloning biases compared to the Sanger method.

Errors are not same with the two high throughput sequencing methods.



Aury et al. (2008) *BMC Genomics*


Resequencing applications

Goal: to analyze various genomes compared to a reference one.

Search for polymorphisms and structural variants in populations, mutation identification in biotechnology, organism evolution analyses, cell differentiation along time, ancient DNA discovery...


Metagenomics: genome characterisation in samples.

A wide range of applications: characterise pathogen micro-organisms in patient tissues, definition of the species found in environment samples, understand species evolution...




From JGI DOE

Functional genomic applications



Wold et al. (2008) *Nat. Methods*

Perspectives



Shendure & Aiden (2012) *Nat. Biotech.*

From bulk to localised single-cell

